


Mediflex Announces expansion of Table Mounted Retractor Line

August 1, 2010. Islandia New York.

Mediflex Surgical Products has announced a series of new products that will be part of their successful **Universal Round Bar Table Mounted Retractor Line**. The new items allow a surgeon to set up a *Bilateral System* by using 2 Field Posts mounted on opposite sides of an OR Table. A new Universal Bilateral Crossbar Frame connects the Field Posts and 2 new Lateral Arms can be connected to either the Field Posts or to the Crossbar. Also added were new Post Couplers that will attach the 1/2" diameter Crossbar Frame and the Lateral Arms to the 1" diameter Field Post.



"We felt it important to offer this version of our Universal Round Bar System particularly for our International Markets," said Dan Adler, International Sales Manager at Mediflex. "This type of set-up is widely used in Asia and elsewhere for upper abdominal procedures such as Liver Transplant," he added. "Additional information is available from our Local Distributors or from www.mediflex.com."

Mediflex is a Division of Flexbar Machine Corp and has specialized in designing, developing and manufacturing unique surgical devices for over 40 years. Their product lines include Bookler® Table Mounted Retractor Systems, Mini-Bookler® Table Mounted Retractor Systems, Davida™ Neurosurgical Retractor Systems, Flex Arm™ Surgical Holders and Positioners, StrongArm™ Surgical Holders and Positioners and the Cutting Edge™ System of Disposable Laparoscopic Instrument Tips.