



Mediflex Announces The Mediflex **Elevator™** Adjustable Height Field Post for Table Mounted Retractor Systems

October 5, 2012. Islandia, NY

Jon Adler, President of Mediflex Surgical Products, announced the availability of the Mediflex **Elevator** Adjustable Height Field Post for use with Table Mounted Retractor Systems. The announcement was made earlier this week during the American College of Surgeons annual meeting in Chicago, IL.

"The Field Post has remained virtually unchanged in the years since John Bookwalter and David Adler collaborated on the development of the Bookwalter System over 30 years ago." said Mr. Adler. "The **Elevator** Field Post allows the surgeon to raise or lower the height of the post to what is needed for a particular procedure. And, this can minimize vertical interference with members of the surgical team."

The device is available in 2 versions. Ref. 72011-EL can be set in 1.25" (3.2cm) increments from 9" to 14" (23cm - 35cm.) Ref. 72009-EL goes from 12" to 18" (30cm - 45cm) in the same increments. With their 1" (2.54cm) diameter, the posts can be adapted to any Bookwalter™ Styled System using a 1" Post Coupler. It can also be used with Mediflex's **Universal Round Bar System** which has also used a Bookwalter Styled Post.

In addition to its height adjustability, the devices are about 30% lighter than the original 1" diameter Field Posts.

Hospitals and Surgi-Centers should contact Mediflex or their Mediflex Distributor for additional information or to schedule an evaluation.

Mediflex is a Division of Flexbar Machine Corp and has specialized in designing, developing and manufacturing unique surgical devices for over 40 years. Their product lines are marketed world-wide and include Bookler®, Universal Round Bar, Mini-Bookler and Davida™ Table Mounted Retractor Systems; Flex Arm™, StrongArm™ and StrongArm **PLUS™** Surgical Holders and Positioners and a wide range of reusable and disposable laparoscopic instruments and devices. The vast majority of Mediflex Products are manufactured in the USA.

For additional information or to receive a free catalog contact:

Mediflex Surgical Products
250 Gibbs Road
Islandia, NY 11749
1-631-582-6424
Email: sales@mediflex.com
Website: www.mediflex.com